

Diagnostic solutions for human infectious diseases

# ELISA kits

for antibody detection of infectious diseases

With more than 20 years of experience, Vircell offers ELISA single tests for detecting specific antibodies against Herpes virus, TORCH and a wide range of infectious agents causing respiratory diseases and sexually transmitted diseases, among others, in human serum/plasma samples.

## Product features

- No sample pre-dilution required.
- Samples and controls are equally processed to compensate pipetting variability.
- Suitable for automated ELISA systems.
- Colour-coded plates with individual break-apart wells.
- Coloured, ready-to-use liquid reagents.
- High performance and guaranteed stability with lyophilized conjugates when necessary.
- Additional semi quantification protocols where quantification is relevant for diagnosis.
- IgM detection in indirect ELISA: absorption of IgG directly in the well using VIRCELL ELISA SORBENT.
- G/M presentations allowing the performance of all possible combinations of IgG and IgM determinations.



# ELISA kits

for antibody detection of infectious diseases

## Single Tests

Product description	Cat. No.	Pack size
ADENOVIRUS ELISA IgG/IgM	G/M1004	96 tests
BRUCELLA ELISA IgG	G1003	96 tests
BRUCELLA ELISA IgM	M1006	96 tests
CHAGAS ELISA IgG+IgM	T1020	96 tests
CHIKUNGUNYA ELISA IgG*	G1010	96 tests
CHIKUNGUNYA ELISA IgM*	M1010	96 tests
CHLAMYDIA ELISA IgG/IgM	G/M1003	96 tests
CHLAMYDIA TRACHOMATIS ELISA IgA	A1017	96 tests
CHLAMYDIA TRACHOMATIS ELISA IgG/IgM	G/M1017	96 tests
CHLAMYDOPHILA PNEUMONIAE ELISA IgA	A1007	96 tests
CHLAMYDOPHILA PNEUMONIAE ELISA IgG (includes calibrators)	G1007	96 tests
CHLAMYDOPHILA PNEUMONIAE ELISA IgM	M1007	96 tests
COXIELLA BURNETII ELISA IgG	G1001	96 tests
COXIELLA BURNETII ELISA IgM	M1001	96 tests
CYTOMEGALOVIRUS ELISA IgG (includes calibrators)	G1004	96 tests
CYTOMEGALOVIRUS ELISA IgM CAPTURE	M1004	96 tests
DENGUE ELISA IgG	G1018	96 tests
DENGUE ELISA IgM CAPTURE	M1018	96 tests
EPSTEIN-BARR EA ELISA IgG	G1205	96 tests
EPSTEIN-BARR EBNA ELISA IgG	G1105	96 tests
EPSTEIN-BARR VCA ELISA IgG	G1005	96 tests
EPSTEIN-BARR VCA ELISA IgM	M1005	96 tests
HELICOBACTER PYLORI ELISA IgA (includes calibrators)	A1022	96 tests
HELICOBACTER PYLORI ELISA IgG (includes calibrators)	G1022	96 tests
HERPES SIMPLEX 1 ELISA IgG/IgM (gG1 recombinant)	G/M1012	96 tests
HERPES SIMPLEX 2 ELISA IgG/IgM (gG2 purified)	G/M1013	96 tests
HERPES SIMPLEX 1+2 ELISA IgG/IgM	G/M1016	96 tests
HYDATIDOSIS ELISA IgG	G1006	96 tests
INFLUENZA A ELISA IgG/IgM	G/M1007	96 tests
INFLUENZA B ELISA IgG/IgM	G/M1008	96 tests
LEGIONELLA PNEUMOPHILA SG 1 ELISA IgG	G1000	96 tests
LEGIONELLA PNEUMOPHILA SG 1 ELISA IgM	M1000	96 tests
LEGIONELLA PNEUMOPHILA SG 1-6 ELISA IgG+IgM	T1001	96 tests
LEISHMANIA ELISA IgG+IgM	T1000	96 tests
LEPTOSPIRA ELISA IgM*	M1003	96 tests
MEASLES ELISA IgG/IgM	G/M1001	96 tests
MUMPS ELISA IgG/IgM	G/M1014	96 tests
MYCOPLASMA PNEUMONIAE ELISA IgG	G1002	96 tests
MYCOPLASMA PNEUMONIAE ELISA IgG/IgM	G/M1015	96 tests
MYCOPLASMA PNEUMONIAE ELISA IgM	M1002	96 tests

Product description	Cat. No.	Pack size
PARAINFLUENZA 1 ELISA IgG/IgM	G/M1009	96 tests
PARAINFLUENZA 2 ELISA IgG/IgM	G/M1010	96 tests
PARAINFLUENZA 3 ELISA IgG/IgM	G/M1011	96 tests
PARVOVIRUS ELISA IgG	G1031	96 tests
PARVOVIRUS ELISA IgM CAPTURE	M1031	96 tests
RESPIRATORY SYNCYTIAL VIRUS ELISA IgG/IgM	G/M1006	96 tests
RICKETTSIA CONORII ELISA IgG/IgM	G/M1005	96 tests
RUBELLA ELISA IgG (includes calibrators)	G1026	96 tests
RUBELLA ELISA IgM CAPTURE	M1026	96 tests
SYPHILIS ELISA IgG	G1060	96 tests
SYPHILIS ELISA IgG+IgM	T1060	96 tests
TETANUS ELISA IgG (includes calibrators)	G1008	96 tests
TOXOPLASMA ELISA IgG (includes calibrators)	G1027	96 tests
TOXOPLASMA ELISA IgM CAPTURE	M1027	96 tests
VARICELLA-ZOSTER ELISA IgG/IgM	G/M1002	96 tests
ZIKA ELISA IgG	G1023	96 tests
ZIKA ELISA IgM	M1023	96 tests

\* Coming soon

## Multitests for antibody detection

Enzyme-linked immunosorbent assays for the simultaneous screening of different ethiological agents in human serum/plasma. The ELISA multitests contain control sera and conjugates for the different types of antigens.

Product description	Cat. No.	Pack size
PNEUMOBACT ELISA IgG <i>Legionella pneumophila</i> sg 1, <i>Mycoplasma pneumoniae</i> , <i>Coxiella burnetii</i> and <i>Chlamydomphila pneumoniae</i>	G1040	96 tests
PNEUMOBACT ELISA IgM <i>Legionella pneumophila</i> sg 1, <i>Mycoplasma pneumoniae</i> , <i>Coxiella burnetii</i> and <i>Chlamydomphila pneumoniae</i>	M1040	96 tests
PNEUMOVIRUS ELISA IgG RSV, Influenza A, Influenza B, Parainfluenza 1, Parainfluenza 2 and Parainfluenza 3	G1050	96 tests

## Other ELISA products

Product description	Cat. No.	Pack size
ELISA SORBENT (2,4 ml)	S001	92 tests
ELISA SORBENT (100 ml)	S001-100	4000 tests
WASH BUFFER FOR ELISA KITS (500 ml 20x)	B0001	3000 tests
SERUM DILUENT FOR ELISA KITS (500 ml)	B0003	5000 tests

All Vircell ELISA IgA and IgM references are adjusted for performance with Vircell ELISA SORBENT only.