

ELISA KITS

For antibody detection of infectious diseases

With more than 20 years of experience, Vircell offers ELISA kits for detecting specific antibodies against Herpes virus, TORCH and a wide range of infectious agents causing respiratory diseases and sexually transmitted diseases, among others, in human serum/plasma samples.

MAIN BENEFITS

- Ready-to-use reagents.
- No sample predilution required.
- Colour-coded plates with individual break-apart wells.
- Semiquantitative protocols available for some parameters.
- μ -Capture technique for TORCH, Dengue and Parvovirus IgM determinations.
- G/M kits allow for customized testing of IgG and/or IgM.

A SOLUTION FOR ALL

- Bacterial Respiratory Infections
- Viral Respiratory Infections
- TORCH Infections
- Vaccine-preventable Infections
- Vector-borne Infections
- Sexually Transmitted Infections
- Zoonotic Infections

OTHER THINGS THAT MATTER

- Suitable for all automated ELISA systems.
- In-house production of high-quality antigens.
- Featured in more than 120 scientific publications.
- Vircell takes part in External Quality Control Programs.



ELISA KITS FOR A WIDE RANGE OF INFECTIOUS DISEASES

BACTERIAL RESPIRATORY INFECTIONS

G/M1003	CHLAMYDIA ELISA IgG/IgM
A1007	CHLAMYDOPHILA PNEUMONIAE ELISA IgA
G1007	CHLAMYDOPHILA PNEUMONIAE ELISA IgG (incl. calibrators)
M1007	CHLAMYDOPHILA PNEUMONIAE ELISA IgM
G1001	COXIELLA BURNETII ELISA IgG
M1001	COXIELLA BURNETII ELISA IgM
G1000	LEGIONELLA PNEUMOPHILA SG 1 ELISA IgG
M1000	LEGIONELLA PNEUMOPHILA SG 1 ELISA IgM
T1001	LEGIONELLA PNEUMOPHILA SG 1-6 ELISA IgG+IgM
G1002	MYCOPLASMA PNEUMONIAE ELISA IgG
G/M1015	MYCOPLASMA PNEUMONIAE ELISA IgG/IgM
M1002	MYCOPLASMA PNEUMONIAE ELISA IgM
M1040	PNEUMOBACT ELISA IgM <i>L. pneumophila</i> sg 1, <i>M. pneumoniae</i> , <i>C. burnetii</i> and <i>C. pneumoniae</i>

VIRAL RESPIRATORY INFECTIONS

G/M1004	ADENOVIRUS ELISA IgG/IgM
G/M1007	INFLUENZA A ELISA IgG/IgM
G/M1008	INFLUENZA B ELISA IgG/IgM
G/M1009	PARAINFLUENZA 1 ELISA IgG/IgM
G/M1010	PARAINFLUENZA 2 ELISA IgG/IgM
G/M1011	PARAINFLUENZA 3 ELISA IgG/IgM
G/M1006	RESPIRATORY SYNCYTIAL VIRUS ELISA IgG/IgM

TORCH INFECTIONS

G1004	CYTOMEGALOVIRUS ELISA IgG (includes calibrators)
M1004	CYTOMEGALOVIRUS ELISA IgM CAPTURE
G/M1012	HERPES SIMPLEX 1 ELISA IgG/IgM (gG1 recombinant)
G/M1016	HERPES SIMPLEX 1+2 ELISA IgG/IgM
G/M1013	HERPES SIMPLEX 2 ELISA IgG/IgM (gG2 purified)
G1026	RUBELLA ELISA IgG (includes calibrators)
M1026	RUBELLA ELISA IgM CAPTURE
G1027	TOXOPLASMA ELISA IgG (includes calibrators)
M1027	TOXOPLASMA ELISA IgM CAPTURE

VACCINE-PREVENTABLE INFECTIONS

G/M1001	MEASLES ELISA IgG/IgM
G/M1014	MUMPS ELISA IgG/IgM
G1008	TETANUS ELISA IgG (includes calibrators)
G/M1002	VARICELLA-ZOSTER ELISA IgG/IgM

VECTOR-BORNE INFECTIONS

T1020	CHAGAS ELISA IgG+IgM
G1010	CHIKUNGUNYA ELISA IgG
M1010	CHIKUNGUNYA ELISA IgM
G1018	DENGUE ELISA IgG
M1018	DENGUE ELISA IgM CAPTURE
T1000	LEISHMANIA ELISA IgG+IgM
G/M1005	RICKETTSIA CONORII ELISA IgG/IgM
G1011	ROSS RIVER VIRUS ELISA IgG
M1011	ROSS RIVER VIRUS ELISA IgM
G1023	ZIKA ELISA IgG
M1023	ZIKA ELISA IgM

SEXUALLY TRANSMITTED INFECTIONS

A1017	CHLAMYDIA TRACHOMATIS ELISA IgA
G/M1017	CHLAMYDIA TRACHOMATIS ELISA IgG/IgM
G1060	SYPHILIS ELISA IgG
T1060	SYPHILIS ELISA IgG+IgM

ZOONOTIC INFECTIONS

G1003	BRUCELLA ELISA IgG
M1006	BRUCELLA ELISA IgM
G1006	HYDATIDOSIS ELISA IgG

OTHER INFECTIONS

G1205	EPSTEIN-BARR EA ELISA IgG
G1105	EPSTEIN-BARR EBNA ELISA IgG
G1005	EPSTEIN-BARR VCA ELISA IgG
M1005	EPSTEIN-BARR VCA ELISA IgM
A1022	HELICOBACTER PYLORI ELISA IgA (incl. calibrators)
G1022	HELICOBACTER PYLORI ELISA IgG (incl. calibrators)
M1003	LEPTOSPIRA ELISA IgM
G1031	PARVOVIRUS ELISA IgG
M1031	PARVOVIRUS ELISA IgM CAPTURE

OTHER ELISA PRODUCTS

S001	ELISA SORBENT (2 x 1.2 ml)
S001-100	ELISA SORBENT (100 ml)
B003	SERUM DILUENT FOR ELISA KITS (500 ml)
B001	WASH BUFFER FOR ELISA KITS (500 ml 20x)

Vircell, S.L.

Parque Tecnológico de la Salud, Avicena 8. 18016 Granada, Spain
+34 958 441 264 | info@vircell.com | www.vircell.com